
Tulip GEMS

TULIP SOFT HARVEST SYSTEM

INSTRUCTIONS FOR USE

ENGLISH

Indications for Use – Tulip Soft Harvest System:

The Tulip Soft Harvest System consists of single-use instrumentation/devices intended for closed-system harvesting and processing of lipoaspirate in medical procedures involving the harvesting, concentrating, and transferring of autologous adipose tissue. The devices are intended for use in the following surgical specialties when the transfer of adipose tissue is desired: orthopedic surgery, arthroscopic surgery, neurosurgery, gastrointestinal and affiliated organ surgery, dermatological surgery and plastic and reconstructive surgery when aesthetic body contouring is intended. Only legally marketed accessory items, such as syringes, should be used with the system. If harvested fat is to be transferred, the harvested fat is only to be used without any additional manipulation.

CAUTION:


- Tulip devices are intended for use with Becton Dickinson (BD) luer lock syringes.
- Verify kit contents.
- Before opening a sterile product, visually inspect the seal around the pouch for any openings.
- Do not use if the pouch is compromised or damaged.
- Open product in a sterile field in accordance with appropriate office, clinical or operating room procedures.
- Use of these devices is limited to those physicians who have attained proficiency in general surgical procedures including suction lipoplasty.

WARNINGS

SINGLE-USE DEVICES ARE NOT TO BE REPROCESSED OR REUSED

- Devices should be used with extreme caution in patients with chronic medical conditions such as heart, lung, or circulatory disease, diabetes, or obesity.
- The volume of blood loss and endogenous body fluid loss may adversely affect intra and/ or postoperative hemodynamic stability. The capability of providing adequate, timely fluid replacement is essential for patient safety.
- When finished using the device, dispose of properly in accordance with biohazard regulations.

Tulip® Single-Use Convenience Kit:

Single-use Cannulas and Accessories are packaged together in various combinations to create Tulip® Single-Use Convenience Kits. Refer to Kit label or included guide to determine the combination of instruments included. The item lot [LOT] number, REF [REF] number, expiration date , and description is located on each individual pouch as well as the case or kit label.

Note: Where referenced, a companion illustration guide is available to enhance the written instructions found in this IFU. These pictorial guides are found on Tulip's website at www.tulipaesthetics.com and may also be included as an insert with that item or kit.

Tulip® Single-Use Cannula:

Single-use cannulas are available in a variety of tip styles, diameters and lengths. They are supplied in pre-sterilized, Tyvek pouches. The item lot [LOT] number, REF [REF] number, expiration date , and description are located on its pouch. The cannula family consists of the following devices:

- Infiltrators
- Injectors
- Spoon Injectors
- Carraway Harvesters
- Miller Harvesters
- Sforza Harvesters
- Sorensen Harvesters
- Tonnard Harvesters

Intended Use – Cannulas:

Used for the Infiltration, Harvesting and Injection of autologous adipose tissue and/or fluids.

Instructions for Use – Cannulas:


Prep: Verify the correct style, diameter and length for the desired procedure. Open by peeling appropriate end of the pouch to remove cannula.

- Tightly screw the cannula onto a sterile BD luer lock syringe. (Syringe size is determined by user)
- Prime the Cannula:
 - Infiltration cannulas – prime the cannula with tumescent solution to eliminate air space.
 - Harvesting cannulas – prime the cannula with saline to eliminate air space and increase suction efficacy, then "Instructions for Use – Snaplocks" section.
 - Injection cannulas – prime the injectors with injectable adipose to eliminate air space

Troubleshooting – Cannulas:

- Remnants of liquids and oils on the mating surfaces of syringes and cannulas can affect the seal and contribute to loss of suction. If loss of suction occurs, 1) ensure connection to syringe is not loose (tighten if required). If suction loss continues, 2) replace with new, dry cannula.
- Note:** Removing/reattaching syringe(s) from a cannula more than 3 times in a single procedure may negatively affect the seal.
- If cannula becomes clogged, insert the cannula tip into a splash-proof container and clear the contents from the cannula with quick, short, in-and-out strokes. Flush with saline.

Tulip® Single-Use Accessories:

Single-use Accessories are supplied in pre-sterilized, Tyvek pouches. The item lot [LOT] number, REF [REF] number, expiration date , and description are located on its pouch. The accessory family consists of the following devices:

- Anaerobic Transfers
- NanoFat Transfers/Cartridges
- SnapLocks
- Syringe Caps

Intended Use – Accessories:

All accessories are used for aid in harvesting, cleaning, and sizing of adipose.

Intended Use – Anaerobic Transfers:

Used for the transfer and sizing of autologous adipose tissue and/or fluid. Available in three sizes:

- 2.4 mm – Used for transfer/emulsification of adipose (injectable by 19 GA needles or larger)
- 1.4 mm – Used for sizing adipose (injectable by 21 GA needles or larger)
- 1.2 mm – Used for sizing adipose (injectable by 23 GA needles or larger)

Note: see the *Injectability Guide* for details

Instructions for Use – Anaerobic Transfers:

Prep: Remove the device from its sterile packaging.

- Attach BD luer lock syringes to each open end of the device.
- Test for a secure connection, making sure the device is seated tightly onto each syringe.
- Transfer adipose/fluid from one syringe to another.
- Repeat transfer step #4 30 times.

Note: *NanoTransfer Set Generation 2 pictorial companion guide* available on www.tulipaesthetics.com.

Intended Use – NanoFat Transfers:

Used for the sizing of autologous adipose tissue to a level at which it can pass through a 30-gauge hypodermic needle.

Note: requires pre-sized adipose using the Anaerobic Transfers. [See Intended Use – Anaerobic Transfers and Anaerobic Transfers Instructions for Use]

Instructions for Use – NanoFat Transfers:

Prep: Remove the device from its sterile packaging. Remove two BD luer lock syringes from their sterile packaging. (Syringe size is determined by user)

- Fill one syringe with pre-sized adipose; empty the other of air completely.
- Attach the empty syringe to the output port of the device; attach the syringe filled with adipose to the input port of the device
- Test for a secure connection, making sure each syringe is tightly sealed to the device.
- Place the device onto a flat surface.
- Transfer adipose from input syringe to the output syringe one time.

Note: *NanoTransfer Set Generation 2 illustration companion guide* available on www.tulipmedical.com.

Intended Use – NanoFat Cartridges:

Used for the sizing of autologous adipose tissue to a level at which it can pass through a 30-gauge hypodermic needle. Used with the NFT housing.

Note: requires pre-sized adipose using the Anaerobic Transfers - NanoTransfer Set Generation 2 pictorial companion guide available on www.tulipaesthetics.com.

Instructions for Use – NanoFat Cartridges:

Prep: Remove the device from its sterile packaging.

- Place cartridge into sterile housing.
 - Note:** The housing is a non-sterile device that requires sterilization before use [see IFU for reusable devices] (*TRUPI 2017 000*).
- Place o-ring onto cartridge.
- Close housing and secure.
- Remove two BD luer lock syringes from their sterile packaging; fill one syringe with pre-sized adipose; empty the other of air completely.
- Attach the empty syringe to the output port of the device. Attach the syringe filled with pre-sized adipose the input port of the device.
- Test for a secure connection, making sure the housing is seated tightly onto each syringe.
- Transfer adipose from the input syringe to the output syringe one time.

Note: *NanoTransfer Set Generation 1 pictorial companion guide* available on www.tulipaesthetics.com.

Intended Use – SnapLocks:

Used for locking/holding BD syringes (under vacuum) at specified volume levels. Fit 20 mL, 50 mL & 60 mL BD luer lock syringes.

Instructions for Use – SnapLocks:


Prep: Remove device from its sterile packaging. Remove appropriate BD luer lock syringe from its sterile packaging.

- Pull the syringe plunger out to its full extent. Orient the snap lock so the plunger clip is adjacent to the base of the plunger.
- Place the Johnnie or Miller Snap onto the plunger head. Using firm pressure, press the snap lock until it clips onto the plunger head.
- The Johnnie or Miller Snap is now securely attached to the plunger and ready to be inserted into the syringe barrel.
- Push the plunger/lock back into the syringe, using your finger to push down each tooth as needed to clear the syringe barrel.
- Once the plunger is seated, the syringe is ready for use.
- Attach the desired suction cannula to the syringe. Prime the cannula with saline to eliminate the air space. And insert the cannula into the soft tissue. (Use small puncture openings rather than incisions for optimal closed system vacuum conditions.)
- Grip the plunger and pull the plunger back until teeth click and anchor the plunger in the desired position.

Note: a vacuum will not be created in the syringe if the cannula tip is not inside body. **Note:** see the *Tulip Single-Use Syringe Lock pictorial companion guide* available on www.tulipaesthetics.com.

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