

Ultrasound

The SonicTouch Series


SonicTouch™ 3 Frequency Ultrasound



State-of-the-Art

ULTRASOUND WITH SMARTHEADS

Tulips' new SmartHeads automatically calibrate to any 50 SeriesPLUS device. SmartHeads in 2, 5, and 10cm² sizes, offer the added benefits of head warming and coupling indication.



3-FREQUENCY ULTRASOUND


The only ultrasound devices in the industry offering 1, 2, and 3 MHz frequencies for the greatest flexibility in depth of treatment.

SonicTouch™ 850 Combination Unit **plus**



3- AND 6-CHANNEL MODEL

The only 3- and 6-channel devices available. All channels allow fully-independent treatment setups. Treatment options include Interferential, premodulated, biphasic, high volt Russian and microcurrent, e-stim and three frequency ultrasound.



The patented TARGET feature remains the standard for focusing the treatment where it's needed. Simply glide your finger across the TARGET touch pad to move the center of interference to the sites of your patient's pain.

TWO-YEAR WARRANTY

Tulip backs these reliable, high quality devices with a full two-year warranty.

TULIP MEDICAL

TO ORDER CALL 888-558-8547

www.tulipmedical.com