



Adipose Substrate Downsizing Using Tulip Anaerobic & Sizing Transfers

Adipose Tissue



Anaerobic Transfer
2.4mm



2.4mm Microfat
Injectable Through
19 Gauge Cannula
or Needle



Sizing Transfer
1.4mm



1.4mm Microfat
Injectable Through
20 Gauge Cannula
or Needle



Sizing Transfer
1.2mm



1.2mm Microfat
Injectable Through
23 Gauge Cannula
or Needle



Nano Transfer



Nanofat
Injectable Through
30 Gauge Needle




For more instructions visit tulipmedical.com/getstarted

MM to Gauge Conversion Guide

- 0.3mm–25GA
- 0.5mm–23GA
- 0.7mm–20GA
- 0.9mm–19GA
- 1.2mm–18GA
- 1.4mm–17GA
- 1.6mm–16GA
- 2.1mm–14GA
- 2.4mm–13GA
- 3.0mm–11GA
- 3.7mm–09GA
- 4.6mm–07GA

Add Fat to Support Natural Healing



23OR003R1