




Harvest Cannula Guide

The Tulip single-use small-volume procedure kit comes with everything you need for an office-based aesthetic procedure and comes with your choice of one of four different harvesters. Choosing the right harvester for your procedure will ensure you achieve the best results for your patients.



Carraway

Miller

Sorensen

Tonnard

	Carraway	Miller	Sorensen	Tonnard
Office-Based Procedures	5	4	3	1
Patient Comfort (Non-bruising)	5	5	4	2
Harvest Speed	3	5	4	4
Use in Fibrous Fat	3	3	4	5
Ease of Use	5	4	3	3
Bloodless Harvest	5	5	4	2

Use our guide below to find the best harvester for your practice. Ratings are based on how you intend to use the harvester. Ranked on a scale of 1 – 5 with 5 being the best for the intended use.



Scan the QR code to Reorder Single-Use Tulip Nano Small Volume Kits