


Recommended Cleaning Products

We recommend using these products to care for your reusable instruments:



Tulip Autoclave Tray

Available in Sizes: Medium, Large

Scan to Purchase



Tulip SoftPicks

Available in Sizes: Micro, Small, Medium

Sold in Packs of 10

Scan to Purchase



Tulip CellFriendly Cleaning Brushes

Available in Sizes: Extra Small, Small, Medium, Medium - Long, Large, Extra-Large - Long

Sold in Packs of 10

Scan to Purchase



Recommended Cleaning Solutions:

Miltex® EZ-Zyme® All Purpose Enzyme Cleaner

Miltex® Surgical Instrument Cleaner



WARNING: DO NOT USE ALKALINE-BASED DETERGENTS OR CLEANING AGENTS TO CLEAN TULIP CANNULAS OR ACCESSORIES. USE OF ALKALINE-BASED DETERGENTS VOIDS TULIP MEDICAL PRODUCT WARRANTY

Tulip® Instrument Care Guide

WARNING: DO NOT USE ALKALINE-BASED DETERGENTS OR CLEANING AGENTS TO CLEAN TULIP CANNULAS OR ACCESSORIES. USE OF ALKALINE-BASED DETERGENTS VOIDS TULIP MEDICAL PRODUCT WARRANTY

To maximize the life of your reusable instruments, always follow the Reusable Instrument Care & Reprocessing Guide.

Scan the QR Code to view the
Reusable Instrument Care & Reprocessing Guide



Tulip has validated its reusable cannulas, injectors, and infiltrators for up to 50 autoclave cycles and strongly recommends replacement after 50 autoclave cycles.

Questions? Call Tulip Medical Products at 858.270.5900 or email sales@tulipmedical.com