

Tulip[®]

Product Catalog



23MCPTP1001R3

Delivering the Benefits of Fat Since 1991
Visit nanofat.com

Table of Contents

Tulip Single-Use Products & Accessories	2-4
Tulip Reusable Products & Accessories	5-9
Syringes & Syringe Accessories	10
Specialty Instruments	11-12
Recommended Sets	13-19
Tulip Hub Styles and the Tulip Difference	20
Contact Information	Back Cover

Tulip Single-Use Products & Accessories



GEMS Tumescant Infiltrator™

This spiral multi-port instrument allows even distribution of tumescent solution.

Single-Use box of 10.

Item Code:	Stem Style:	Size:
DC1SLINF1412	Infiltrator	2.1mm (14G) x 12cm
DC1SLINF1415	Infiltrator	2.1mm (14G) x 15cm
DC1SLINF1420	Infiltrator	2.1mm (14G) x 20cm



GEMS Carraway Harvester™

This spiral three-port design is safe and effective for multi-plane harvesting. The Carraway harvester is commonly used for small-volume procedures performed using local anesthesia.

Single-Use box of 10.

Item Code:	Stem Style:	Size:
DC1SLCAR1315	Carraway	2.4mm (13G) x 15cm
DC1SLCAR1410	Carraway	2.1mm (14G) x 10cm
DC1SLCAR1415	Carraway	2.1mm (14G) x 15cm



GEMS Miller Speed Harvester™

A six-port harvester that harvests maximum tissue volume in minimal time.

Single-Use box of 10.

Item Code:	Stem Style:	Size:
DC1SLMIL1415	Miller Speed	2.1mm (14G) x 15cm



GEMS Sorensen Harvester™

Designed by Dr. Jesper Sorensen, this is a less aggressive multiport microfat harvester with beveled micro ports.

Single-Use box of 10.

Item Code:	Stem Style:	Size:
DC1SLSOR1320	Sorensen	2.4mm (13G) x 20cm

Tulip Single-Use Products & Accessories



GEMS Tonnard Harvester™

Designed by Dr. Patrick Tonnard, this unique multiport harvester, with elevated micro ports, harvests microfat efficiently.

20cm=20-ports

Single-Use box of 10.

Item Code:	Stem Style:	Size:
DC1SLTON1320	Tonnard	2.4mm (13G) x 20cm



GEMS Tulip Injector™

Tulip's side port, blunt tip design passes through soft tissue effortlessly and maintains the integrity of the vascular structure. A workhorse for safe, small-volume injection.

Single-Use box of 10.

Item Code:	Stem Style:	Size:
DC1SLINJ1709	Injector	1.4mm (17G) x 09cm
DC1SLINJ1807	Injector	1.2mm (18G) x 07cm
DC1SLINJ1905	Injector	0.9mm (19G) x 05cm
DC1SLINJ2004	Injector	0.7mm (20G) x 04cm
DC1SLINJ2306	Injector	0.5mm (23G) x 06cm
DC1SLINJ2505	Injector	0.3mm (25G) x 05cm



GEMS Spoon Tip Injector™

A blunt spoon-shaped tip designed with a larger port for rapid injection of larger volumes - your workhorse for breast and buttock injection.

Single-Use box of 10.

Item Code:	Stem Style:	Size:
DC1SLSPO1415	Spoon Tip Injector	2.1mm (14G) x 15cm



Tulip Single-Use Products & Accessories



GEMS Johnnie Snap™ For 60cc and 20cc BD Syringes Single-Use Box of 5

Patented GEMS™ Johnnie Snap 60cc and 20cc syringe locks are designed to maintain vacuum at variable pressure levels, from low to high. With a newly designed plunger clip, this lock becomes one with the plunger delivering effortless vacuum pressure. The Johnnie Snap delivers incremental, low-pressure vacuum with 2 lock positions for 20cc syringes and 3 lock positions for 60cc syringes.

Item Code:

A. **DJS220** Johnnie Snap 20cc

B. **DJS360G2** Johnnie Snap 60cc ThreeSnap



GEMS Miller OneSnap™ For 60cc BD Syringes Single-Use Box of 5

Patented Single-Use Miller OneSnap 60cc syringe locks are designed to maintain a level vacuum. With a newly designed plunger clip, this lock becomes one with the plunger delivering effortless vacuum pressure.

The Miller OneSnap delivers a high-pressure vacuum.

Item Code:

DMS160G2 Miller Snap 60cc OneSnap

Tulip Reusable Products & Accessories



Tumescent Infiltrator™

Also available in single-use

The spiral, multi-port instrument allows even distribution of tumescent solution.

Item Code:	Stem Style:	Hub Style:	Size:
INFSL1320CF	Infiltrator	SuperLuerLok Hub	2.4mm (13G) x 20cm
INFSL1325CF	Infiltrator	SuperLuerLok Hub	2.4mm (13G) x 25cm
INFSL1412CF	Infiltrator	SuperLuerLok Hub	2.1mm (14G) x 12cm
INFSL1415CF	Infiltrator	SuperLuerLok Hub	2.1mm (14G) x 15cm
INFSL1420CF	Infiltrator	SuperLuerLok Hub	2.1mm (14G) x 20cm
INFLL1325CF	Infiltrator	Luer Lock	2.4mm (13G) x 25cm
INFLL1415CF	Infiltrator	Luer Lock	2.1mm (14G) x 15cm
INFLL1420CF	Infiltrator	Luer Lock	2.1mm (14G) x 20cm
INF601125CF	Infiltrator	60cc Toomey Hub	3.0mm (11G) x 25cm



Carraway Harvester™

Also available in single-use

This spiral three-port design is safe and effective for multi-plane harvesting. The Carraway harvester is commonly used for small-volume procedures performed using local anesthesia.

Only available on the SuperLuerLok™ Hub.

Item Code:	Stem Style:	Hub Style:	Size:
CARSL1315CF	Carraway	SuperLuerLok Hub	2.4mm (13G) x 15cm
CARSL1320CF	Carraway	SuperLuerLok Hub	2.4mm (13G) x 20cm
CARSL1410CF	Carraway	SuperLuerLok Hub	2.1mm (14G) x 10cm
CARSL1415CF	Carraway	SuperLuerLok Hub	2.1mm (14G) x 15cm



Flap™

The Flap is a multipurpose dissector, harvester and injector.

Only available on the SuperLuerLok™ Hub.



Item Code:	Stem Style:	Hub Style:	Size:
FLPSL1315CF	Flap	SuperLuerLok Hub	2.4mm (13G) x 15cm
FLPSL1415CF	Flap	SuperLuerLok Hub	2.1mm (14G) x 15cm



Mercedes Harvester™

This three-port design harvests from all directions.

Item Code:	Stem Style:	Hub Style:	Size:
MERSL1320CF	Mercedes	SuperLuerLok Hub	2.4mm (13G) x 20cm
MER600925CF	Mercedes	Toomey Hub	3.7mm (9G) x 25cm
MER601125CF	Mercedes	Toomey Hub	3.0mm (11G) x 25cm

Tulip Reusable Products & Accessories



Miller Speed Harvester™

Also available in single-use

A six-port harvester harvest maximum tissue volume in minimal time.

Item Code:	Stem Style:	Hub Style:	Size:
MILSL1315CF	Miller Speed Harvester	SuperLuerLok Hub	2.4mm (13G) x 15cm
MILSL1320CF	Miller Speed Harvester	SuperLuerLok Hub	2.4mm (13G) x 20cm
MILSL1415CF	Miller Speed Harvester	SuperLuerLok Hub	2.1mm (14G) x 15cm
MIL601125CF	Miller Speed Harvester	60cc Toomey Hub	3.0mm (11G) x 25cm



Sorensen Harvester™

Also available in single-use

Designed by Dr. Jesper Sorensen, this is a less aggressive multiport microfat harvester with beveled micro ports.

Only available on the SuperLuerLok™ Hub.

Item Code:	Stem Style:	Hub Style:	Size:
SORSL1320CF	Sorensen Harvester	SuperLuerLok Hub	2.4mm (13G) x 20cm



Standard™

The single-port design conservatively harvests from one direction.

Item Code:	Stem Style:	Hub Style:	Size:
STDSL1315CF	Standard	SuperLuerLok Hub	2.4mm (13G) x 15cm
STDSL1320CF	Standard	SuperLuerLok Hub	2.4mm (13G) x 20cm
STDSL1415CF	Standard	SuperLuerLok Hub	2.1mm (14G) x 15cm
STD600925CF	Standard	60cc Toomey Hub	3.7mm (9G) x 25cm
STD601125CF	Standard	60cc Toomey Hub	3.0mm (11G) x 25cm



Tonnard Harvester™

Also available in single-use

Designed by Dr. Patrick Tonnard, this unique multiport harvester, with elevated sharp micro ports, harvests microfat efficiently. **15cm and shorter = 12-ports**
20cm and longer = 20-ports

Item Code:	Stem Style:	Hub Style:	Size:
TONSL1120CF	Tonnard Harvester	SuperLuerLok Hub	3.0mm (11G) x 20cm
TONSL1315CF	Tonnard Harvester	SuperLuerLok Hub	2.4mm (13G) x 15cm
TONSL1320CF	Tonnard Harvester	SuperLuerLok Hub	2.4mm (13G) x 20cm
TONSL1415CF	Tonnard Harvester	SuperLuerLok Hub	2.1mm (14G) x 15cm
TONLL1120CF	Tonnard Harvester	Luer Lock	3.0mm (11G) x 20cm
TON601120CF	Tonnard Harvester	60cc Toomey Hub	3.0mm (11G) x 20cm

Tulip Reusable Products & Accessories



Toledo V-Dissector™

Designed by Dr. Luis Toledo, this multipurpose instrument undermines and releases connective tissue and injects new tissue into the space.

Now available portless—see below

Item Code:	Stem Style:	Hub Style:	Size:
VDISL1415CF	V-Dissector	SuperLuerLok Hub	2.1mm (14G) x 15cm
VDISL1420CF	V-Dissector	SuperLuerLok Hub	2.1mm (14G) x 20cm
VDISL1607CF	V-Dissector	SuperLuerLok Hub	1.6mm (16G) x 07cm



Fakih Portless V-Dissector™

The Fakih Portless V-Dissector is a specialized tool that is used to prepare dense and fibrous tissue for infiltration and tissue extraction. The cutting notch is very useful for undermining in facial procedures and in areas where the tissue must be broken down due to acne scarring, stretch marks, and other fibrous deposits.

Item Code: VDINP1609



Bensimon Micro Sculptor

Designed by Dr. Richard Bensimon, this modification of the Tonnard harvester is ideal for facial liposuction and inducing skin contraction. The small holes ensure a smooth contour, and the small diameter enables entry with a needle puncture leaving no mark. This allows for multiple entries for precise contouring of the neck and gynecomastia. The Bensimon Micro Sculptor is a 12-Hole SuperLuerLok CellFriendly cannula that is barbed and extra toothy.

Item Code:	Hub Style:	Size:
BENSL1415CF	SuperLuerLok Hub	2.1mm (14G) x 15cm
BENSL1315CF	SuperLuerLok Hub	2.4mm (13G) x 15cm

Tulip Reusable Products & Accessories



Tulip Injector™

Also available in single-use

Tulip's side port, blunt tip design passes through soft tissue effortlessly and maintains the integrity of the vascular structure. A workhorse for safe, small-volume injection.

Item Code:	Stem Style:	Hub Style:	Size:
INJSL1709CF	Injector	SuperLuerLok Hub	1.4mm (17G) x 09cm
INJSL1807CF	Injector	SuperLuerLok Hub	1.2mm (18G) x 07cm
INJSL1905CF	Injector	SuperLuerLok Hub	0.9mm (19G) x 05cm
INJSL2004CF	Injector	SuperLuerLok Hub	0.7mm (20G) x 04cm
INJSL2005CF	Injector	SuperLuerLok Hub	0.7mm (20G) x 05cm
INJLL1807CF	Injector	Luer Lock	1.2mm (18G) x 07cm
INJLL1905CF	Injector	Luer Lock	0.9mm (19G) x 05cm
INJLL2004CF	Injector	Luer Lock	0.7mm (20G) x 04cm



Spoon Tip Injector™

Also available in single-use

A blunt, spoon-shaped tip designed with a larger port for rapid injection of larger volumes — your workhorse for breast and buttock injection.

Item Code:	Stem Style:	Hub Style:	Size:
SPOSL1120CF	Spoon Tip Injector	SuperLuerLok Hub	3.0mm (11G) x 20cm
SPOSL1315CF	Spoon Tip Injector	SuperLuerLok Hub	2.4mm (13G) x 15cm
SPOSL1320CF	Spoon Tip Injector	SuperLuerLok Hub	2.4mm (13G) x 20cm
SPOSL1415CF	Spoon Tip Injector	SuperLuerLok Hub	2.1mm (14G) x 15cm
SPOLL1320CF	Spoon Tip Injector	Luer Lock Hub	2.4mm (13G) x 20cm
SPO601120CF	Spoon Tip Injector	60cc Toomey Hub	3.0mm (11G) x 20cm



Curved Spoon Tip Injector™

Tulip's uniquely curved injectors give the surgeon precise control of the injection plane. Turned downward, the curved injectors volumetrically fill in the deep plane. Turned up the injectors glide against the superficial plane to stimulate regenerative results.

Item Code:	Stem Style:	Hub Style:	Size:
CSPSL1415CF	Curved Spoon Injector	SuperLuerLok Hub	2.1mm (14G) x 15cm

Tulip Reusable Products & Accessories

CellFriendly Instrument Cleaning Brushes™

These brushes are specially designed to clean CellFriendly Cannulas. Soft proprietary nylon bristles clean the lumen without damage to the CellFriendly surface.

Only use Tulip brushes and SoftPicks on Tulip cannulas to maintain the CellFriendly benefits.

Sold in packages of 10.

Brush XS – Fits 0.9mm and 1.2mm. Length is 15cm

Brush SM – Fits 1.4mm and 1.6mm. Length is 20cm

Brush Med – Fits 2.1mm, 2.4mm, 3.0mm. Length is 30cm

Brush Med/Long – Fits 2.1mm, 2.4mm, 3.0mm handles. Length is 60cm

Brush Large – Fits 3.7mm and 4.6mm. Length is 40cm

Brush Long XL – Fits 3.7mm and 4.6mm for handles. Length is 60cm


Tulip SoftPicks™

Cleaning even the smallest reusable CellFriendly micro injectors is quick, easy, and gentle. Guaranteed not to score, abrade, or contaminate. Tulip Cannula Cleaning SoftPicks are the ideal single-use cleaning solution for smaller-diameter cannulas and microinjectors. Autoclavable for use in the sterile field, if desired.

Sold in packages of 10.

- **Micro** – For cleaning 0.7mm & 0.9mm cannulas
Item Code: TSP 19/20
- **Small** – For cleaning 1.2mm, 1.4mm & 1.6mm cannulas
Item Code: TSP 16/17/18
- **Medium** – For cleaning 2.1mm, 2.4mm & 3.0mm cannulas
Item Code: TSP 11/13/14

Choosing the Right Tools to Clean Cannulas or Injectors



DISCLAIMER: CANNULA/INJECTOR DIAMETERS ARE BASED ON TULIP CANNULA/INJECTORS ONLY
NOTE: IF CANNULA OR INJECTOR HAS HANDLE ATTACHED PLEASE KEEP FULL HANDLE AND CANNULA/INJECTOR LENGTH IN MIND


23ER007R1

Syringes and Syringe Accessories



Tulip Syringes

- A. 60cc Toomey Monoject Syringe - Box of 20
Item Code: SYR60BX
- B. 60cc Luer BD Syringe - Box of 40
Item Code: SYR60BDBX
- C. 20cc Luer BD Syringe - Box of 48
Item Code: SYR20BDBX
- D. 10cc Luer BD Syringe - Box of 100
Item Code: SYR10BDBX
- E. 1cc Luer BD Polycarbonate Syringe - Box of 100
Item Code: SYR1LLBDBX



Syringe Tip Caps

Keep the contents of your syringe sealed.

Options for syringe caps available:

- A. Luer Tip Caps
Sterile, single-use 5 packs of 10 caps each (50 single-use caps total)
Item code: CAPSTD-05
- B. Luer Tip caps for Centrifuge Sterile, Single-Use
5 packs of 10 caps each (50 single-use caps total)
Item code: CAPSLCENT-05
- C. Toomey Tip Caps Non-Sterile, Reusable
For 60cc Toomey Tip syringes (1 pack of 5 caps)
Item code: CAP60-05



SnapLok™ - Reusable

These syringe SnapLoks fit into the plunger flanges of the liposuction syringe. When the plunger is pulled out the lock snaps open on the barrel lip holding the liposuction syringe under vacuum for CellFriendly fat extraction.

Available in the following sizes:

- Item Code:**
- SNAP10** 10cc (for 10cc BD syringes)
- SNAP20** 20cc (for 20cc BD syringes)
- SNAP60** 60cc (for 60cc Monoject syringes)



Johnnie Lok™

The reusable Johnnie Lok holds the plunger under vacuum pressure for easy syringe harvesting. The plunger is drawn back to the desired pressure level then twisted and locked securely.

Available in the following sizes:

- Item Code:**
- LOK10BD** 10cc BD (for 10cc BD syringes)
- LOK20BD** 20cc BD (for 20cc BD syringes)
- LOK60BD** 60cc BD (for 60cc BD syringes)
- LOK60** 60cc Monoject (for 60cc Monoject syringes)

Specialty Instruments

NanoSpin™

Advanced small-volume tissue processing

- Designed as a closed system, works with our Tulip Crowns and syringe caps (purchased separately)
- Speed up your workflow and increase efficiency
- Decrease contamination
- Preprogrammed settings for proper separations of adipose, PRP, and PRF
- Autoclavable Syringe Sleeves
- Rotor holds eight 10 mL syringes
- 500-3500 rpm speed range
- 1,534g Max G-force
- 12V with a DC motor
- Made in the USA



Item Code: NanoSpin G2

You may also purchase the 10cc Tulip Crowns (See below)

Tulip Crown™

U.S. Patent Number 10039886

The Crown's patented design keeps the syringe contents in a virtually anaerobic atmosphere during centrifuging. Simply attach a Crown to the syringe plunger, harvest, then twist the plunger off. The Crown also has a bonus, self-sealing extractor port for cleanly removing contents from any layer of centrifuged tissue.



Item Code:

CROWN10BD

CROWN20BD

CROWN60

Available in sizes:

10cc BD *(compatible with NanoSpin)

20cc BD

60cc Monoject

Tulip Gold Standard Autoclave Tray™

Tulip's Gold Standard autoclave trays are custom-designed to protect, house, and store Tulip reusable sets and individual Tulip instruments. (Pictured is Medium Tray)



Available in the following sizes:

Medium Tray

Item Code: TrayMed

For Micro and Facial Cannulas

Dimensions: 15.24 x 25.40 x 4.45 cm

Upper Mat for Medium Tray

Code: UpperMatTrayGDMed

Large Tray *Limited Supply*

Item Code: TrayGSLG

For Facial and Body Cannulas

Dimensions: 15.24 x 40.64 x 4.45 cm

Specialty Instruments



Tulip Power Adapter™

The Tulip Power Adapter modifies the syringe barrel for use with a liposuction aspirating machine.

Size:
20cc BD
60cc Monoject

Item Code:
PA20BD
PA60



Pocar™

The Pocar is a percutaneous punch that is used to make a circumferential incision prior to infiltration and tissue extraction. It is designed to maintain a tight seal while minimizing scarring.

Size:
Large 4.6mm
Small 3.0mm

Item Code:
POCAR4.6
POCAR3.0



Tulip Z Stand™

The ideal syringe stand for keeping your decanted tissue organized and at the ready.

The autoclavable Tulip Z Stand is available in 2 sizes:

Size:
Large | fits 20cc–60cc syringes
Small | fits 1cc–12cc syringes

Item Code:
SYRSTDLG
SYRSTDSM



Decanting Stand™

Built-in anaerobic transfers are used to prepare and transfer extracted tissue for re-injection. The decant stand design also allows for a natural gravity decant of extracted tissue.

Sizes
60cc Toomey to Luer Lock
Luer Lock to Luer Lock

Item Code:
AT60LLDS
ATLLDS

Recommended Sets

Tissue Processing Microfat and Nanofat

Tulip's Nanofat System is designed around Tulip's clinically validated protocol for producing nanofat. Tulip designed this protocol and the instrumentation in collaboration with the pioneers of nanofat, Drs. Patrick Tonnard and Alexis Verpaele of Gent, Belgium.

Characteristics of Nanofat:

1. Nanofat must be consistently injectable with a 30G needle
2. Nanofat must remain viable within its native vascular structure and concentrates



Tulip Nanofat Set™ - Reusable

Included in the set:

- 1 each Tulip NanoTransfer™ Generation I
- 2 boxes of 5 Sterile NanoTransfer Cartridges
- 1 each Reusable Anaerobic Transfer 2.4mm
- 1 each Reusable Sizing Transfer 1.4mm
- 1 each Reusable Sizing Transfer 1.2mm

Item Code: TGSNANO-DNTC



Tulip Nanofat Kit™ - Single-Use

This single-use patent-pending closed system helps size adipose tissue so that it is injectable with 27G and 30G needles.

Kit for 1 Patient Includes:

- (Kit for 5 Patients Available with amounts of 5 each per item)
- 1 each Single-Use NanoTransfer™ Generation 2
- 1 each GEMS™ Anaerobic Transfer 2.4mm
- 1 each GEMS™ Sizing Transfer 1.2mm
- 1 each GEMS™ Sizing Transfer 1.4mm

Item Code:

DNTS-1P for 1 Patient

DNTS-5P for 5 Patient



NanoTransfer™ Generation I Cartridges

Pre-sterilized single-use cartridges for use with the Generation I NanoTransfer

Sold in boxes of 5

Item Code: DNTC-05

Recommended Sets

Tissue Processing Microfat and Nanofat




The Micro Kit™ - Single-Use

GEMS™ Single-use Anaerobic Sizing Transfers are used to size down adipose tissue into desired aliquot sizes of Microfat. They also play an important role in Tulip's proprietary Nanofat System™. Transfers can be used as a standalone system for sizing down adipose or can be used in conjunction with the Tulip NanoTransfer (Gen 1 or Gen 2) to create NanoFat. Each of the GEMS™ Transfers is also available for purchase individually in boxes of 10 (Codes in parentheses next to each item).

Kit Includes:

- 5 each GEMS™ Anaerobic Sizing Transfer 2.4mm (DATLLLL2.4-10)
- 5 each GEMS™ Sizing transfer 1.4mm (DATLLLL1.4-10)
- 5 each GEMS™ Sizing transfer 1.2mm (DATLLLL1.2-10)

Item Code: DMICRO-05



The Micro Set™ - Reusable

Reusable Anaerobic Sizing Transfers are used to size down adipose tissue into desired aliquot sizes of Microfat. They also play an important role in Tulip's proprietary Nanofat System™. Transfers can be used as a standalone system for sizing down adipose or can be used in conjunction with the Tulip NanoTransfer (Gen 1 or Gen 2) to create NanoFat. Each Reusable Transfer is also available for purchase individually (code is in parentheses next to each item).

Set Includes:

- 1 each Reusable Anaerobic Transfer 2.4mm (ATLLLL2.4)
- 1 each Reusable Sizing Transfer 1.4mm (ATLLLL1.4)
- 1 each Reusable Sizing Transfer 1.2mm (ATLLLL1.2)

Item Code: MicroSet

Available in Single-Use (as seen above)

Recommended Sets

Tulip Nano Small Volume Kit Single-Use

Tulip Nano Small Volume Kit™ – Single-Use

The reusable Tulip Nano Small Volume Kit is your workhorse for facial sculpting and rejuvenation using microfat and nanofat. It is also ideal for other small-volume procedures such as hand rejuvenation, scar revision, wound healing, décolletage rejuvenation, and restoration of the fat pads in the feet.



Tulip Nano Small Volume Kit™ Single-Use

- 1 each Tumescant Infiltrator SuperLuerLok 2.1mm x 20cm or 2.1 mm x 15cm (Depends on Harvester Size)
 - 2 each Harvester of your choice:
 - Carraway Harvester SuperLuerLok 2.1mm x 15cm
 - Miller Harvester SuperLuerLok 2.1mm x 15cm
 - Tonnard Harvester SuperLuerLok 2.4mm x 20cm
 - Sorensen Harvester SuperLuerLok 2.4mm x 20cm
 - 1 each Tulip Injector SuperLuerLok 0.7mm x 4cm
 - 1 each Tulip Injector SuperLuerLok 0.9mm x 5cm
 - 1 each Tulip Injector SuperLuerLok 1.2mm x 7cm
 - 1 each Miller OneSnap 60cc or Snap 20 (BD)
 - 1 each NanoTransfer Gen 2
 - 2 each Tulip Anaerobic Transfers 2.4mm
 - 1 each Tulip Sizing Transfers 1.4mm
 - 1 each Tulip Sizing Transfers 1.2mm
 - 2 each Harvesting Syringes BD 60cc syringes
 - 2 each Harvesting Syringes BD 10cc syringes
 - 4 each Injection Syringes BD 1cc syringes
- (Packaged as a single procedure kit)



Available in a variety of configurations:

Harvester Options:

Harvesters:	Item Code:
GEMS™ Sorensen Harvester	DPK-SOR-TNF
GEMS™ MillerSpeed Harvester	DPK-MIL-TNF
GEMS™ Tonnard Harvester	DPK-TON-TNF
GEMS™ Carraway Harvester	DPK-CAR-TNF

Recommended Sets

Tulip Nano Small Volume Kit Reusable

Tulip Nano Small Volume Set™ – Reusable

The reusable Tulip Nano Small Volume Set is your workhorse for facial sculpting and rejuvenation using microfat and nanofat. It is also ideal for other small-volume procedures such as hand rejuvenation, scar revision, wound healing, décolletage rejuvenation, and restoration of the fat pads in the feet.



Set Includes:

- 1 each Tumescant Infiltrator SuperLuerLok 2.1mm x 20cm or 2.1mm x 15cm (Depends on Harvester Choice)
- 1 each Harvester of your choice:
 - Carraway Harvester SuperLuerLok 2.1mm x 15cm
 - Miller Harvester SuperLuerLok 2.1mm x 15cm
 - Sorensen Harvester SuperLuerLok 2.4mm x 20cm
 - Tonnard Harvester SuperLuerLok 2.4mm x 20cm
- 1 each Tulip Injector SuperLuerLok 0.7mm x 4cm
- 1 each Tulip Injector SuperLuerLok 0.9mm x 5cm
- 1 each Tulip Injector SuperLuerLok 1.2mm x 7cm
- 1 each SnapLok 60cc
- 1 each NanoTransfer
- 2 each Anaerobic Transfer Luer to Luer 2.4mm
- 1 each Tulip Sizing Transfers 1.4mm
- 1 each Tulip Sizing Transfers 1.2mm
- 2 boxes of 5 Sterile NanoTransfer Cartridges (Packaged as procedure set)



Available in a variety of configurations

Harvester Options:

Harvesters:

Carraway Harvester
Miller Speed Harvester
Sorensen Harvester
Tonnard Harvester

Item Code:

TGS-CAR-TNF
TGS-MIL-TNF
TGS-SOR-TNF
TGS-TON-TNF

Recommended Sets

Tulip Nano SoftBreast Kits



Tulip Nano SoftBreast™ Kit™ Single-Use

This single-use procedure kit for personalized Breast Enhancement, with or without implants, using microfat and nanofat for aesthetic and reconstructive applications.

Item Code: DPK-SFB-TNF

- 1 each GEMS™ Tumescant Infiltrator SuperLuerLok 2.1mm x 15cm
- 2 each GEMS™ Miller Harvester SuperLuerLok 2.1mm x 15cm
- 1 each GEMS™ Spoon Tip Injector SuperLuerLok 2.1mm x 15cm
- 1 each GEMS™ Injector SuperLuerLok 1.4mm x 9cm
- 1 each GEMS™ Injector SuperLuerLok 1.2mm x 7cm
- 2 each GEMS™ Miller OneSnap 60cc
- 1 each NanoTransfer Generation 2
- 2 each GEMS™ Anaerobic Transfers 2.4mm
- 1 each GEMS™ Sizing Transfers 1.4mm
- 1 each GEMS™ Sizing Transfers 1.2mm
- 4 each Harvesting Syringes BD 60cc
- 6 each Harvesting Syringes BD 10cc
- 6 each Injection Syringes BD 1cc
- (Packaged as procedure set)



Tulip Nano SoftBreast™ Set™ Reusable

This single-use procedure set for personalized Breast Enhancement, with or without implants, using microfat and nanofat for aesthetic and reconstructive applications.

Item Code: TGSSB-TNF

- 1 each Tumescant Infiltrator SuperLuerLok 2.1mm x 15cm
- 1 each Miller Speed Harvester 2.1mm x 15cm
- 1 each Spoon Tip Injector 2.1mm x 15cm
- 1 each Tulip Injector 1.4mm x 9cm
- 1 each Tulip Injector 1.2mm x 7cm
- 1 each SnapLok 60cc
- 1 each NanoTransfer
- 2 boxes of 5 Sterile NanoTransfer Cartridges
- 2 each Reusable Anaerobic Transfer 2.4mm
- 1 each Reusable Sizing Transfer 1.4mm
- 1 each Reusable Sizing Transfer 1.2mm
- (Packaged as procedure set)

Recommended Sets

Doctor Sets



Marten Clinic Set™

The Marten Set was designed and inspired by Dr. Timothy Marten of San Francisco, CA. This set for the discriminating fat grafting facelift surgeon is a comprehensive collection of CellFriendly instrumentation and accessories for refined facial sculpting using microfat and nanofat.

Item Code: MAR-TNF

- 1 each Tumescant Infiltrator SuperLuerLok 2.1mm x 20cm
- 1 each Carraway Harvester SuperLuerLok 2.1mm x 15cm
- 1 each Carraway Harvester SuperLuerLok 2.4mm x 20cm
- 1 each Tonnard Harvester SuperLuerLok 2.1mm x 15cm
- 2 each Tulip Injector SuperLuerLok 0.7mm x 4cm
- 2 each Tulip Injector SuperLuerLok 0.7mm x 5cm
- 2 each Tulip Injector SuperLuerLok 0.9mm x 5cm
- 2 each Tulip Injector SuperLuerLok 1.2mm x 7cm
- 1 each V-Dissector SuperLuerLok 1.6mm x 7cm
- 1 each NanoTransfer
- 2 boxes of 5 Sterile NanoTransfer Cartridges
- 2 each Reusable Anaerobic Transfer 2.4mm
- 1 each Reusable Sizing Transfers 1.4mm
- 1 each Reusable Sizing Transfers 1.2mm
- (Packaged as procedure set)



Tonnard/Verpaele Set™


From the Pioneers of Nanofat, this reusable and CellFriendly set, designed and recommended by Dr. Patrick Tonnard and Dr. Alexis Verpaele, is perfect for fine facial sculpting and contouring.

Item Code: TVS-TNF

- 1 each Tumescant Infiltrator SuperLuerLok 2.1mm x 20cm
- 1 each Tonnard Harvester SuperLuerLok 2.4mm x 20cm
- 1 each Tulip Injector SuperLuerLok 0.7mm x 5cm
- 1 each Tulip Injector SuperLuerLok 0.9mm x 5cm
- 1 each SnapLok 60cc
- 1 each SnapLok 20cc
- 1 each NanoTransfer
- 2 boxes of 5 Sterile NanoTransfer Cartridges
- 2 each Reusable Anaerobic Transfer 2.4mm
- 1 each Reusable Sizing Transfers 1.4mm
- 1 each Reusable Sizing Transfers 1.2mm
- (Packaged as procedure set)

Recommended Sets

Doctor Sets



Ghavami Gluteal Set™

This CellFriendly set, inspired and recommended by Dr. Ashkan Ghavami of Beverly Hills, CA., is designed for expert body contouring. This versatile set is ideal for Mommy Makeovers, breast enhancement, gluteal augmentation, reshaping the mid-torso, and producing a unified, curvaceous shape from neck to knees.

Item Code: GAV-TNF

- 1 each Mercedes Harvester 60cc 3.7mm x 20cm
 - 1 each Mercedes Harvester 60cc 3.7mm x 25cm
 - 1 each SpoonTip Injector 60cc 3.0mm x 15cm
 - 1 each SpoonTip Injector 60cc 3.0mm x 20cm
 - 1 each NanoTransfer
 - 2 boxes of 5 Sterile NanoTransfer Cartridges
 - 1 each Reusable Anaerobic Transfer 2.4mm
 - 1 each Reusable Sizing Transfers 1.4mm
 - 1 each Reusable Sizing Transfers 1.2mm
 - 1 each Decant Stand Anaerobic Transfer 60cc to Luer
- (Packaged as procedure set)



Nabil Fakh Master Face Set™

Curated by Dr. Nabil Fakh, this set is designed for facial lipofilling using Microfat, SNIF (Sharp Needle Intradermal Fat Grafting), SNIE (Sharp Needle Intradermal Emulsified Fat Grafting), and Nanofat. This all-in-one set equips you to achieve a natural and harmonious result when sculpting and contouring the face, rejuvenating the skin, and revising scars.

Item Code: FAK

- 1 each Tumescent Infiltrator SuperLuerLok 2.1mm x 20cm
 - 1 each Tonnard Harvester SuperLuerLok 2.4mm x 20cm
 - 1 each Tulip Injector SuperLuerLok 0.7mm x 4cm
 - 1 each Tulip Injector SuperLuerLok 0.9mm x 5cm
 - 1 each Tulip Injector SuperLuerLok 1.2mm x 7cm
 - 1 each Fakh Portless V-Dissector 1.6mm x 09cm
 - 1 each SnapLok 20cc
 - 1 each JohnnieLok BD 10cc
 - 1 each NanoTransfer
 - 2 boxes of 5 Sterile NanoTransfer Cartridges
 - 1 each Reusable Anaerobic Transfer 2.4mm
 - 1 each Reusable Sizing Transfers 1.4mm
 - 1 each Reusable Sizing Transfers 1.2mm
 - 1 each Syringe Stand Large
 - 1 each Anaerobic Transfer Luer to Luer Decant Stand
- (Packaged as procedure set)

Hub Information & The Tulip CellFriendly™ Difference



60cc Hub™

Tulip SuperLuerLok™

Luer Lock


Standard Instrument Dimensions

GAUGE MICRO

25GA	0.3mm 0.51mm OD(0.31mm ID)
23GA	0.5mm 0.64mm OD(0.43mm ID)
20GA	0.7mm 0.91mm OD(0.66mm ID)
19GA	0.9mm 1.07mm OD(0.81mm ID)
18GA	1.2mm 1.27mm OD(0.98mm ID)
17GA	1.4mm 1.47mm OD(1.21mm ID)
16GA	1.6mm 1.65mm OD(1.36mm ID)

MACRO

14GA	2.1mm 2.11mm OD(1.70mm ID)
13GA	2.4mm 2.41mm OD(1.96mm ID)
11GA	3.0mm 3.05mm OD(2.39mm ID)
09GA	3.7mm 3.76mm OD(3.00mm ID)
07GA	4.6mm 4.57mm OD(3.81mm ID)



Get the Tulip Advantage using CellFriendly™ Instruments!

All of Tulip's cannulas and injectors are finished with Tulip's exclusive CellFriendly™ coating. This finish gives the instruments a smooth, lubricious action. Reduced friction creates less stress and damage to the tissue for a better outcome and less fatigue to the surgeon during the procedure.

SUPERCARGE your microfat grafts

Surgeons worldwide are now supplementing their microfat with nanofat to supercharge their injections. These new micro injectors are sized right for injecting nanofat - supercharged microfat. Whether you are using sized-down microfat alone or nanofat - supercharged microfat, these micro injectors are ideal for periorbital, perioral, nose, and forehead applications. They are also excellent for restoring soft tissue in the foot fat pad as well as for female wellness restoration.

Contact Information



Call us at +1 (858) 270-5900 or
Toll-free +1 (800) 325-6526



Visit TulipAesthetics.com to order online or
email us at sales@tulipmedical.com



For international orders, contact
distributor@tulipmedical.com or visit the
website to view the list of available regions
at Tulipmedical.com/distributors

TULIP® and Tulip Nano® are registered trademarks of Black Tie Medical Inc., DBA Tulip Medical Products. © 2023 Black Tie Medical Inc., DBA Tulip Medical Products. All rights reserved.