



nanoTM SPIN

Compact Design

- Great for Office-based Procedures
- Preprogrammed for Adipose Separation
- Process Eight 10 ml syringes (max 80ml)

Price: \$1,700.00



nanoTM SPIN PRO

Maximum Capacity

- Ideal for Frequent Use
- Preprogrammed for Adipose Separation
- Process Six 10ml or 20ml Syringes (max 120ml)
- CE Certified

Price: \$2,500.00

	NANOSPIN TM	NANOSPIN TM PRO
Use Tulip Crowns TM and Syringe Caps to Make it a Closed System	✓	✓
Optimal Separations of Adipose	✓	✓
CE Certified		✓
Made in the USA	✓	✓
Autoclavable Syringe Sleeves	✓	✓
Capacity	8 10ml syringes	6 10ml or 20ml syringes
Speed Range	500-3500 RPM	Max 3900 RPM
Max G-Force	1,534 G	1,850 G
Power	110V/220V	110V/220V
LED Lid Lights for Cycle Status		✓
Locking Lid to Prevent Improper Opening	Manual	Automatic
Interface	Pre-programmed	Pre-programmed
Cycle Memory	✓	✓

Advanced Closed-System Tissue Processing

The Patented Tulip Crown™ May be Small, But it Makes a Big Difference.

- Allows for easy removal of the syringe plunger
- Autoclavable for 50 Autoclave Cycles
- Compatible with both the NanoSpin™ and NanoSpin Pro™
- Prevent unwanted messes
- Decrease contamination
- Use Tulip Crowns with syringe caps for virtually anaerobic tissue processing
- Patented Design fits 10ml or 20ml BD Syringes

Price: \$160.00/each

Crown Sizes Available:

10ml BD (compatible with NanoSpin & Nanospin Pro)

20ml BD (compatible with Nanospin Pro)

